

Total RNA extraction from *Oryza sativa* leaf and *Solanum tuberosum* leaf using Ribospin™ Seed/Fruit

Experimental Conditions

Materials Required

- Ribospin™ Seed/Fruit (50 preps: 317-150)
- TissueLyser II (85300, supplier Q) or another bead beating device
- Liquid nitrogen (LN₂)
- Absolute ethanol (C₂H₆O, CAS No.: 64-17-5, ≥99.0%)
- β-mercaptoethanol (C₂H₆OS CAS No.: 60-24-2, ≥99.0%)
- 1.5 ml or 2.0 ml microcentrifuge tube
- Vortex mixer
- Centrifuge (Max. speed 14,000 rpm or ≥10,000 x g)
- Pipette & sterile pipette tips
- Suitable protector (e.g., lab coat, disposable gloves, goggles, etc.)
- Ice

Sample Information

Sample type:



Oryza sativa leaf



Solanum tuberosum leaf

- Sampling: after collecting fresh leaves, put them in a bag and seal it.
- How to store: store in -70 °C deep freezer
- Homogenizing: TissueLyser II (85300, supplier Q)
- Conditions
 - Sample amount
 - Oryza sativa* leaf: 100 mg
 - Solanum tuberosum* leaf: 80 mg
 - Elution volume: 50 µl

Protocol

Before experiment

1. Before using for the first time, add absolute ethanol (ACS grade or better) into Buffer RBW and RNW as indicated on the bottle.
2. Prepare DNase I reaction mixture just before step 6.
 - Prepare aliquot DNase I and thaw on ice.
 - Mix 2 µl DNase I with 70 µl Buffer DRB.

Sample preparation

1. After measuring 100 mg of *Oryza sativa* leaf sample and 80 mg of *Solanum tuberosum* leaf sample, put them into a 1.5 ml microcentrifuge tube. Then, it is rapidly frozen using LN₂.
2. Grind the frozen sample using TissueLyser II.
 - *Oryza sativa* leaf samples: 30 Hz, 40 sec
 - *Solanum tuberosum* leaf samples: 30 Hz, 20 sec
3. The next step is according to protocol II (not protocol I) of Ribospin™ Seed/Fruit.

Ribospin™ Seed/Fruit Protocol (Protocol II)

* For more details and methods, please refer to [the handbook of Ribospin™ Seed/Fruit](#).

1. Add 500 µl Buffer SL and 5 µl β-mercaptoethanol to the sample and vortex vigorously for 15 sec.
2. Incubate the mixture for 3 min at room temperature and centrifuge the lysate at 13,000 rpm for 1 min. Transfer 300 µl of the supernatant to a new 1.5 ml microcentrifuge tube (not provided).
3. Add 300 µl Buffer ML to the supernatant and vortex vigorously for 15 sec and transfer all of the mixture to EzPure™ Filter (yellow).
4. Centrifuge at 13,000 rpm for 1 min and transfer 500 µl of the pass-through to a new 1.5 ml microcentrifuge tube (not provided).
5. Add 250 µl absolute ethanol to the supernatant and mix it well by inversion. Apply all of the mixture into Column Type F (blue ring) and centrifuge at 13,000 rpm for 1 min.
6. Add 500 µl Buffer RBW to Column Type F and centrifuge at 13,000 rpm for 30 sec.
7. Apply 70 µl DNase I reaction mixture onto the center of Column Type F for gDNA digestion. Incubate for 10 min at room temperature.
8. Add 500 µl Buffer RBW to Column Type F and centrifuge at 13,000 rpm for 30 sec. Add 500 µl Buffer RNW to Column Type F and centrifuge at 13,000 rpm for 30 sec.
9. Centrifuge at maximum speed for an additional 1 min to remove residual wash buffer. Transfer Column Type F to a new 1.5 ml microcentrifuge tube (provided).
10. Add 50 µl Nuclease-free Water to the center of the membrane in Column Type F and centrifuge at 13,000 rpm for 1 min.

Total RNA extraction from *Oryza sativa* leaf and *Solanum tuberosum* leaf using Ribospin™ Seed/Fruit

Result

#	Sample	Conc. (ng/μl)	A _{260/280}	A _{260/230}	Yield (μg)
1	<i>Oryza sativa</i> leaf (100 mg)	183.4	2.13	2.35	9.17
2		173.4	2.14	2.38	8.67
3		164.3	2.15	2.35	8.22
4	<i>Solanum tuberosum</i> leaf (80 mg)	111.9	2.15	2.20	5.56
5		106.9	2.16	2.18	5.35

Table 1. The concentrations, yield and purity of RNA extracted from 100 mg of *Oryza sativa* leaf and 80 mg of *Solanum tuberosum* leaf samples.

※ Absorbance measurement instrument: NanoDrop™ 2000/2000c (ND-2000, supplier T)




Figure 1. The result of electrophoresis of total RNA from 100 mg of *Oryza sativa* leaf and 80 mg of *Solanum tuberosum* leaf samples.

Lane M: GENESTA™ 1 kb DNA Ladder with 5X loading dye (GA-100, GeneAll®, 1 μl loading)

Lanes 1-3: 100 mg of *Oryza sativa* leaf samples (3 μl loading)

Lanes 4-5: 80 mg of *Solanum tuberosum* leaf samples (5 μl loading)

※ Electrophoresis conditions: 1.0% agarose gel (150 V, 17 min)