

# Comparative Evaluation of Genomic DNA Extraction Performance Between Exgene™ Clinic SV mini and Other Commercial DNA Extraction Kits

## Experimental Conditions

### Materials Required

- Exgene™ Clinic SV mini (Cat. No. 108-101 or 108-152)
- 1.5 ml & 2.0 ml microcentrifuge tube
- 50 ml conical tube
- EDTA tube
- Sterilized swab
- Centrifuge & microcentrifuge ( $\leq 15,000 \times g$ )
- Vortex mixer
- Absolute ethanol ( $\geq 99.0\%$ ,  $C_2H_5OH$ , CAS No. 64-17-5)
- 1X PBS (Phosphate Buffered Saline), pH 7.4
- Pipette & sterilized pipette tips
- Suitable protector (e.g., lab coat, disposable gloves, goggles, etc.)

### Sample Information

- Extraction conditions

Sample	Sample Amount	Elution Volume
Human whole blood	200 $\mu$ l	50 $\mu$ l
Saliva (mouthwash)	10 ml	
Buccal swab	1 stick	
Cultured cell (K562)	$5 \times 10^6$ cells	
Urine	1 ml	

### Sample Preparation

#### • Human whole blood

1. Prepare human whole blood in EDTA tube or blood collection tube with anticoagulants mixture.
2. Transfer 200  $\mu$ l of human whole blood to the 1.5 ml microcentrifuge tube and follow the [A. protocol for blood and body fluid/cultured cells using microcentrifuge](#) (Page 18).

#### • Saliva (mouthwash)

1. Collect 10 ml of mouthwash into a 50 ml conical tube and add 5 ml of 1X PBS.
2. Vortex the mixture to mix thoroughly and centrifuge at  $2,000 \times g$  for 5 min. Carefully decant the supernatant.
3. Resuspend the pellets completely in 200  $\mu$ l of 1X PBS and follow the [C. protocol for saliva and mouthwash](#) (Page 24).

#### • Buccal swab

1. Collect the oral epithelial cells using sterile swab.
2. Place the swab into a 2.0 ml microcentrifuge tube and cut off the handle of swab with a sterile sharp blade or cutter.
3. Add 400  $\mu$ l of 1X PBS to the tube and follow the [B. protocol for buccal swab](#) (Page 22).

#### • Cultured cell (K562)

1. Transfer the harvested cells into a 1.5 ml microcentrifuge tube and centrifuge at  $14,000 \times g$  for 1 min.
2. Discard the supernatant and resuspend the cell pellet with 200  $\mu$ l of 1X PBS. Follow the [A. protocol for blood and body fluid/cultured cells using microcentrifuge](#) (Page 18).

#### • Urine

1. Transfer 1 ml of urine to a 1.5 ml microcentrifuge tube and centrifuge at  $6,000 \times g$  above for 2 min.
2. Discard the supernatant and add 200  $\mu$ l of 1X PBS then vortex the tube for 5 sec.
3. Centrifuge at  $6,000 \times g$  above for 2 min. Then discard the supernatant.
4. Follow the [A. protocol for blood and body fluid/cultured cells using microcentrifuge](#) (Page 18).

\* For more details and methods, please refer to the [handbook of Exgene™ Blood/Clinic/Cell SV mini protocol](#).


## Result

Sample	Kit	Yield ( $\mu$ g)	$A_{260/280}$
Human whole blood	Exgene™ Clinic SV mini	<b>5.15</b>	1.86
	Supplier A	<b>2.62</b>	1.90
	Supplier B	<b>2.93</b>	1.83
Saliva (mouthwash)	Exgene™ Clinic SV mini	<b>0.78</b>	2.00
	Supplier A	<b>0.61</b>	2.16
	Supplier B	<b>0.91</b>	2.02
Buccal swab	Exgene™ Clinic SV mini	<b>0.21</b>	2.39
	Supplier A	<b>0.16</b>	2.56
	Supplier B	<b>0.17</b>	3.39
Cultured cell (K562)	Exgene™ Clinic SV mini	<b>2.17</b>	1.89
	Supplier A	<b>1.60</b>	1.90
	Supplier B	<b>3.11</b>	2.00
Urine	Exgene™ Clinic SV mini	<b>0.05</b>	1.21
	Supplier A	<b>0.09</b>	2.17
	Supplier B	<b>0.12</b>	2.14

**Table 1. Absorbance analysis of genomic DNA extracted from different types of samples with Exgene™ Clinic SV mini and other commercial DNA extraction kits**  
Genomic DNA was extracted three times from each of five different samples using Exgene™ Clinic SV mini and genomic DNA extraction kits from Supplier A and Supplier B.

The concentration and purity of the eluate were determined using a spectrophotometer (NanoDrop™ 2000, Supplier T), and the yield was calculated based on the measured values. Each value represents the average of triplicate measurements.






# Comparative Evaluation of Genomic DNA Extraction Performance Between Exgene™ Clinic SV mini and Other Commercial DNA Extraction Kits



**Figure 1. Electrophoresis analysis of genomic DNA extracted using Exgene™ Clinic SV mini and other commercial DNA extraction kits from various sample types**  
 Genomic DNA was extracted from each sample using Exgene™ Clinic SV mini and two other commercial DNA extraction kits. Extracted DNA in triplicate was analyzed by electrophoresis in a 1.2% agarose gel and visualized under UV light.

**• Lane Information**

- M : GENESTA™ 1 kb DNA Ladder (GA-100)
- Lane 1, 6, 11 : DNA eluate from human whole blood
- Lane 2, 7, 12 : DNA eluate from saliva (mouthwash)
- Lane 3, 8, 13 : DNA eluate from buccal swab
- Lane 4, 9, 14 : DNA eluate from cultured cell (K562)
- Lane 5, 10, 15 : DNA eluate from urine



**Figure 2. C<sub>q</sub> value on a real-time PCR (qPCR) amplification curve obtained from DNA extracted from each sample**

Each amplification curve shows real-time PCR (qPCR) amplification of DNA template extracted from each sample using Exgene™ Clinic SV mini and two other commercial DNA extraction kits. The each C<sub>q</sub> value in the table are represented by different colors.

**• PCR Information**

- PCR primer : Human GAPDH
- Real-time PCR (qPCR) instrument : CFX96™ (1855201, Supplier B)
- Real-time PCR (qPCR) kit : RealAmp™ 2X qPCR Master Mix (801-050)