

Comparative Evaluation of Genomic DNA Extraction Performance Between Exgene™ Clinic SV mini and Other Commercial DNA Extraction Kits

Experimental Conditions

Materials Required

- Exgene™ Clinic SV mini (100 preps: 108-101 / 250 preps: 108-152)
- 1.5 ml & 2.0 ml microcentrifuge tube
- 50 ml conical tube
- EDTA tube
- Sterilized swab
- Centrifuge & microcentrifuge ($\leq 15,000 \times g$)
- Vortex mixer
- Absolute ethanol ($\geq 99.0\%$, C_2H_5OH , CAS No.: 64-17-5)
- 1X PBS (Phosphate-buffered saline), pH 7.4 (SM-P04-100)
- Pipette & sterilized pipette tips
- Suitable protector (e.g., lab coat, disposable gloves, goggles, etc.)

Sample Information

- Extraction conditions

Sample	Sample Amount	Elutoin Volume
Human whole blood	200 μ l	50 μ l
Saliva (mouthwash)	10 ml	
Buccal swab	1 stick	
Cultured cell (K562)	5 x 10 ⁶ cells	
Urine	1 ml	

Sample Preparation

* For more details and methods, please refer to [the handbook of Exgene™ Clinic SV mini, midi, maxi](#).

Human whole blood

- Prepare human whole blood in EDTA tube or blood collection tube with anticoagulants mixture.
- Transfer 200 μ l of human whole blood to the 1.5 ml microcentrifuge tube and follow the **A. Protocol for whole blood/body fluid/cultured cells using microcentrifuge (page 10)**.

Saliva (mouthwash)

- Collect 10 ml of mouthwash into a 50 ml conical tube and add 5 ml of 1X PBS.
- Vortex the mixture to mix thoroughly and centrifuge at 2,000 $\times g$ for 5 min. Carefully decant the supernatant.
- Resuspend the pellets completely in 200 μ l of 1X PBS and follow the **C. Protocol for saliva/mouthwash (page 15)**.

Buccal swab

- Collect the oral epithelial cells using sterile swab.
- Place the swab into a 2.0 ml microcentrifuge tube and cut off the handle of swab with a sterile sharp blade or cutter. Add 400 μ l of 1X PBS to the tube and follow the **B. Protocol for buccal swab (page 13)**.

Cultured cell (K562)

- Transfer the harvested cells into a 1.5 ml microcentrifuge tube and centrifuge at 14,000 $\times g$ for 1 min.
- Discard the supernatant and resuspend the cell pellet with 200 μ l of 1X PBS. Follow the **A. Protocol for whole blood/body fluid/cultured cells using microcentrifuge (page 10)**.

Urine

- Transfer 1 ml of urine to a 1.5 ml microcentrifuge tube and centrifuge at 6,000 $\times g$ above for 2 min.
- Discard the supernatant and add 200 μ l of 1X PBS then vortex the tube for 5 sec.
- Centrifuge at 6,000 $\times g$ above for 2 min. Then discard the supernatant.
- Follow the **A. Protocol for whole blood/body fluid/cultured cells using microcentrifuge (page 10)**.

Result

Sample	Kit	Yield (μ g)	A _{260/280}
Human whole blood	Exgene™ Clinic SV mini	5.15	1.86
	Supplier A	2.62	1.90
	Supplier B	2.93	1.83
Saliva (mouthwash)	Exgene™ Clinic SV mini	0.78	2.00
	Supplier A	0.61	2.16
	Supplier B	0.91	2.02
Buccal swab	Exgene™ Clinic SV mini	0.21	2.39
	Supplier A	0.16	2.56
	Supplier B	0.17	3.39
Cultured cell (K562)	Exgene™ Clinic SV mini	2.17	1.89
	Supplier A	1.60	1.90
	Supplier B	3.11	2.00
Urine	Exgene™ Clinic SV mini	0.05	1.21
	Supplier A	0.09	2.17
	Supplier B	0.12	2.14

Table 1. Absorbance analysis of genomic DNA extracted from different types of samples with Exgene™ Clinic SV mini and other commercial DNA extraction kits
Genomic DNA was extracted three times from each of five different samples using Exgene™ Clinic SV mini and genomic DNA extraction kits from supplier A and supplier B. The concentration and purity of the eluate were determined using a spectrophotometer (NanoDrop™ 2000, supplier T), and the yield was calculated based on the measured values. Each value represents the average of triplicate measurements.

Comparative Evaluation of Genomic DNA Extraction Performance Between Exgene™ Clinic SV mini and Other Commercial DNA Extraction Kits




Figure 1. Electrophoresis analysis of genomic DNA extracted using Exgene™ Clinic SV mini and other commercial DNA extraction kits from various sample types
Genomic DNA was extracted from each sample using Exgene™ Clinic SV mini and two other commercial DNA extraction kits. Extracted DNA in triplicate was analyzed by electrophoresis in a 1.2% agarose gel and visualized under UV light.

• Lane information

- M: GENESTA™ 1 kb DNA Ladder (GA-100)
- Lane 1, 6, 11: DNA eluate from human whole blood
- Lane 2, 7, 12: DNA eluate from saliva (mouthwash)
- Lane 3, 8, 13: DNA eluate from buccal swab
- Lane 4, 9, 14: DNA eluate from cultured cell (K562)
- Lane 5, 10, 15: DNA eluate from urine








Figure 2. C_q value on a real-time PCR (qPCR) amplification curve obtained from DNA extracted from each sample

Each amplification curve shows real-time PCR (qPCR) amplification of DNA template extracted from each sample using Exgene™ Clinic SV mini and two other commercial DNA extraction kits. The each C_q value in the table are represented by different colors.

• PCR information

PCR primer: Human GAPDH
Real-time PCR system: CFX96™ System (1855201, supplier B)
qPCR kit: RealAmp™ 2X qPCR Master Mix (801-050)