


## End of run



Once extraction is completed, take Plate/Tube Rack out of the system and collect genomic DNA from 5<sup>th</sup>/11<sup>th</sup> wells depending on where samples were dispensed. Transfer eluent into 1.5 ml or PCR tube.

Amount of eluted buffer will be around 180 µl and it is ready to use.

## Troubleshooting guide

Problem	Causes	Comments and suggestions
Physical damage of the kit	Upside down during transportation may cause beads to stick with sealing film	Spin down the 'Cartridge' by hand before open it.
	Sealing film is detached and reagent is spilled to other wells due to improper storage temperature	Spin down by hand and measure reagent volume with eyes. If reagent volumes are insufficient, extraction efficiency may decrease. Do not use it and contact customer service immediately.
	Broken 'Cartridge'	Broken 'Cartridge' may lead to unfavorable result. Do not use it and contact customer service immediately.
Inappropriate specimen condition	Specimen condition is not favorable due to inappropriate storage condition (ex. coagulation)	Perform vortex and pipetting of specimen. If sample is still coagulated, add a bit of PBS or distilled water and vortex again.
	Specimen condition is not favorable due to inappropriate storage condition (ex. stored in room temperature)	Increase sample volume up to 300 µl, if extraction efficiency is low.
Abnormal extraction	Too much beads left in Elution buffer	If the total nucleic acid density is within the normal range, proceed with the eluted solution. In the case of low total nucleic acid density, transfer the eluted solution to a 1.5 ml tube and centrifuge before use.
		Refer to 'Inappropriate specimen condition' part of trouble shooting if specimen condition is unfavorable, perform extraction again. If the specimen condition and total nucleic acid density are favorable, proceed with the extracted total nucleic acid. In cases where the specimen condition is favorable but the total nucleic acid density is unfavorable, transfer the eluted solution to 1.5 ml tube and centrifuge before use. If the result remains unfavorable, dilute it with distilled water before use.
	Eluted total nucleic acid should not appear transparent or sticky	

## Warnings and precautions

- Intended for research use only.
- Read and follow the manual before using the product.
- Use extracted nucleic acid as soon as possible, if not, keep it at -70°C for long-term storage.
- Be cautious of contaminants such as microorganisms after opening the product.
- Be sure to wear personal protective equipment such as gloves and goggles when using this product and wash hands after handling specimens and reagents.
- Be mindful of contamination with DNase or RNase during product use.
- Store the product at the specified storage temperature and do not use it past its expiration date.
- Read and follow the IFU for the nucleic acid extraction device (GENTI™ Advanced Automated Nucleic Acid Extraction Equipment) used with this product.
- The reagents in this product contain irritants, do not dispose of them with bleach or acids.
- This product is a single use and should not be reused.

\* Any serious incident involving the device is reported to the relevant competent authority in the country where the manufacturer, user and patient are located.

## Storage conditions

- Temperature : Room temperature (15~25°C)
- Humidity : 20~80%



GeneAll Bldg., 303-7 Dongnam-ro, Songpa-gu, Seoul, 05729, Korea  
E-mail : sales@geneall.com  
Tel : 82-2-407-0096 Fax : 82-2-407-0779



A-1201~A-1204, Hanam Techno Valley U1 Center, 947, Hanam-daero, Hanam-si, Gyeonggi-do, 12982, Korea

## GENEALL BIOTECHNOLOGY CO., LTD.

©2024 GeneAll Biotechnology, All rights reserved

Symbol	Used for	Symbol	Used for
	Batch number		Manufacturer
	Catalogue number		Do not re-use
	Consult instructions for use		Date of manufacture
	Caution		Expiry date
	Temperature limit		

2024.06

Ver. 1.2

GeneAll®

Store at room temperature (15~25°C)  
Expiration date : 18 months after manufacture

GENTI™ Advanced  
Plant Kit (Tube Type/Plate Type)

## Description

The GENTI™ Advanced Plant Kit provides efficient plant DNA/RNA co-extraction, in conjunction with magnetic bead-based GENTI™ Advanced Automatic Extraction Equipment.

Available extraction protocols for DNA and RNA cater to both medium throughput (1 to 16 samples) with Tube Type Cartridge and high throughput (up to 32 samples) with the Plate Type Cartridge.

## Intended use

GENTI™ Advanced Plant Kit is formulated to extract both DNA and RNA from a wide range of plant materials like leaf, seed and bark.

The extracted nucleic acid is suitable for various downstream applications such as PCR, RT-PCR, qPCR, qRT-PCR, NGS and other molecular diagnostic tests.

## Kit contents

Components	Quantity	
	904-048A	904-096A
Number of preparation	48 T	96 T
Type	Tube *	Plate
Reagent pre-filled cartridge	6 pks	6 pks
Magnetic rod cover (6 pcs/pk)	4 pks	2 pks
Buffer EL 150 ml	1 ea	2 ea
Buffer BR 150 ml	1 ea	2 ea
Buffer SQ1 150 ml	1 ea	2 ea
Nuclease-free water 1 ml	1 ea	2 ea

\* The Tube Rack, which is an essential accessory for using the Tube Types, is provided with GENTI™ Advanced Automatic Extraction Equipment.

• Reagent pre-filled cartridge (Plate Type)



• Reagent pre-filled cartridge (Tube Type)



• Disposable magnetic rod cover




• Heating blocks



- Conical shape of the Plate/Tube Type Cartridge, magnetic rod cover and heating block  
- Heating block designed with tight fit with conical materials for fast and efficient heat transfer

## Brief workflow



- Lysis well
- 1<sup>st</sup> washing well
- 2<sup>nd</sup> washing well
- 3<sup>rd</sup> washing + magnetic beads well
- Elution well
- 4<sup>th</sup> washing well

- The Tube Type Cartridge is designed for low to medium throughput, less than 8 samples. It is tailored to adjust the number of preps, contributing to efficient reagent saving.
- The Plate Type Cartridge enables the extraction of 16 samples, utilizing six consecutive wells for each extraction. These six wells are arranged horizontally, with each well housing specific reagents for the extraction process.
- Both kits can be used with same hardware allowing the users to switch between the two methods according to the requirements in sample.
- 1<sup>st</sup> (7<sup>th</sup>) wells contain lysis buffer which disrupts cell membranes and binds target DNA and RNA with magnetic beads.
- Magnetic beads, stored in the 4<sup>th</sup> (10<sup>th</sup>) are moved to the 1<sup>st</sup> (7<sup>th</sup>) well upon extraction initiation.
- 2<sup>nd</sup>, 3<sup>rd</sup>, 4<sup>th</sup>, 6<sup>th</sup> (8<sup>th</sup>, 9<sup>th</sup>, 10<sup>th</sup>, 12<sup>th</sup>) wells contain washing buffer in order to remove unwanted cell component and buffers.
- Elution buffer, housed in the 5<sup>th</sup> (11<sup>th</sup>) well facilitates the dissolution of nucleic acid molecules from the magnetic beads.

## ■ Sample preparation for leaf sample

- Grind sample to a fine powder completely using a mortar and pestle under liquid nitrogen. Place up to 100 mg of ground sample into a 1.5 ml microcentrifuge tube (not provided).  
\* Bark samples can be applied up to 1 g.
- Add 200 µl Buffer EL and 200 µl Buffer BR to the sample and vortex vigorously for 15 sec.  
\*\* For starch-enriched samples, add 300 µl of Buffer EL and 300 µl of Buffer BR.  
\*\*\* For bark samples, add 3 ml of Buffer EL and 3 ml of Buffer BR.
- Incubate the mixture for 3 min at room temperature.
- Centrifuge the lysate at 13,000 rpm (≥10,000 x g) for 5 min at 4°C.
- Transfer 200 µl of the supernatant to 1<sup>st</sup> (7<sup>th</sup>) well.
- Load plate on the tray of GENTi™ Advanced Automatic Extraction Equipment.  
\* Note : Ensure that diagonally cut edge of Tube Rack faces the top left of the heating block and that the Tube Rack is placed evenly.
- Insert magnetic rod cover to the end to strip bracket.  
\* Note : Ensure that magnetic rod cover is in the correct position.

## ■ Sample preparation for seed sample (starch-enriched samples)

- Grind sample to a fine powder completely using a mortar and pestle under liquid nitrogen. Place up to 1 g of ground sample into a 15 ml tube (not provided).
- Add 3 ml Buffer SQ1 and vortex vigorously for 15 sec.
- Continue with step 4 in sample preparation for leaf sample.


## ■ Extraction protocol of GENTi™ Advanced Automatic Extraction Equipment

Protocol	Uses and purpose
Fast Protocol (11' 02")	<ul style="list-style-type: none"> <li>Time-saving, high-speed extraction for PCR-ready nucleic acids</li> </ul>
Normal Protocol (19' 50")	<ul style="list-style-type: none"> <li>Standard procedure of nucleic acid extraction</li> <li>Optimized for nucleic acid extraction from a variety of leaf, root, bulb and bark samples</li> </ul>
High Protocol (28' 26")	<ul style="list-style-type: none"> <li>High quality nucleic acid extraction (seed sample specific protocol)</li> <li>Optimized for nucleic acid extraction from starch-enriched seed samples</li> </ul>

## ■ Progression of normal protocol


Step	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Step 7	Step 8	Step 9
Well	4	1	2	3	4	6	6	5	4
Name	Bead TF	Lysis	Wash 1	Wash 2	Wash 3	Wash 4	Dry	Elution	Reclaim
Mixing	0:00	5:00	0:50	0:50	0:30	0:30	0:00	3:10	0:10
Volume	700	1000	700	700	700	700	100	100	50
Block Tm	Off	35°C	Off	Off	Off	Off	Off	85°C	Off

## ■ Preparation of Tube Type Cartridge




Cut the cartridge as much as necessary and install it in the GENTi™ Advanced Tube Rack.

## ■ Preparation of Plate Type Cartridge



## ■ Equipment run



- Turn on the 'Power switch' located on the right rear of the GENTi™ Advanced Automatic Extraction Equipment.
- Touch the 'RUN' button when the home screen appears.
- Select the 'Self Test' at the File Browser menu and then touch the 'RUN' button to run a self-test.
- After self test completes, select the extraction protocol and then touch the 'RUN' button for the protocol operation.  
(Optimal protocol of the three options, 1) Fast <11 min>, 2) Normal <19 min>, 3) High <28 min>)  
\* Note : If 'Warning' message appears, check the system and touch the 'RUN' button again.

## ■ Precautions for use of equipment



### Self-test

The self-test starts automatically when the equipment turns on.

### Mounting the Plate/Tube Type Cartridge

Gently mount the Tube Rack on the equipment with caution.

### Mounting magnetic rod cover

When installing the magnetic rod cover, push it to the end of magnetic rod cover slot.