



Pharmacidal

A bacteria prevention spray for incubators and work benches in cell culture and molecular biology labs

Cat. No.:	IC-110100	1 liter, spray bottle
	IC-110100-B	100ml, spray bottle
	IC-110100-L	250ml, spray bottle
	IC-110100-G	5 liter, refill

The problem of contamination in incubators and/or sterile work-benches is often a serious one, leading to extensive damage. Pharmacidal solution prevents contamination with and growth of fungi (and spores), bacteria (and spores) and viruses (including HIV and Hepatitis B).

The active ingredients are quarternary benzylammonium compounds, and the solution does not contain mercury, formaldehyde, phenol or alcohol.

Pharmacidal is non-toxic and biodegradable. It is non-corrosive. Pharmacidal was also found it to be fully compatible with common work surfaces.

Recommended use:

- Spray incubators once every 2 weeks.
- Spray sterile benches once a day or preferably before using the work area.
- The surfaces to be disinfected should be completely wet by spraying. Allow to air-dry. The drying time is the reaction time!

Instructions for Use